

Diagnostics Plans



Outline



- # Interface and Integration
- # The Goal and the Obstacles
- # System Plans by Priority
- # ‘Diagnostic Event’ Triggers and Buffering
- # Unified Data Display

Interface and Integration



⌘ Diagnostics are maturing

- Major systems are in place and shaken out
- Effort shifting to Interface and Integration

⌘ Interface to Operations

- Unified Data Displays - Injection, Ramp, Store
- Minimize need for Documentation and Specialists

⌘ Integration During L/P² and Machine Development

- Instruments run from the Sequencer, or
- Instruments integrated with Applications where possible, or
- Instruments operable from MCR w/o specialists

The Goal and the Obstacles



- # How we (Diagnostics) understand our job
- # Goal - to assist in providing L and P^2
- # Obstacles - limits are defined by
 - Control of Tune, Chromaticity, and Coupling
 - Injection Efficiency - matching, damping, IBS, e-cloud/vacuum, beam-beam, dynamic aperture,...
 - Ramp Reliability - IBS, e-cloud/vacuum, beam-beam, instabilities, transition, ...
 - L/P^2 Lifetime at Store - IBS, beam-beam, e-cloud/vacuum, dynamic aperture, resonance compensation, nonlinearities,...
 - Short lifetime requires good injection efficiency and ramp reliability

System Plans by Priority



- ⌘ BLM - recalibrate thresholds
- ⌘ BPM - Todd separate talk this session
- ⌘ WCM/DCCT - no major changes, FFT (split?)
- ⌘ ARTUS - in this talk
- ⌘ PLL - separate talk this session
- ⌘ IPM - Roger separate talk this session; nestled IPM and the optical IPM
- ⌘ Schottky's - in this talk
- ⌘ AC Dipole - in this talk
- ⌘ Dampers - Angelika and Mike separate talks this session
- ⌘ Coherence Monitor - Mike separate talk this session
- ⌘ Head-Tail, Buttons, QMM - in this talk
- ⌘ Stochastic Cooling - good data from last run. Plan?

ARTUS



Version 2 Analog Front End

- Easier and more reliable timing
- Enables another coherence monitor?
- Enables switch to all large β BPMs

All large β BPMs

- Helps compensate for loss of 2m of Kicker to Damper
- Requires additional hardware and effort from Controls

Preamps in the tunnel?

LF and HF Schottky



Summary of Last Run - Schottky did not fulfill its potential

- Saturation - due beam steering and effect of 200MHz RF
- Frequency agility - available local oscillators, mixers, digitizers
- Software (version control and LabVIEW idiosyncrasies)
- GBytes of HF spectra archived - will be examined with refined software

Plans for next run


- Solve aperture problems around IP2 (also important for IPM)
- Position control of Moveable BPM
- Improved frequency agility - NCOs, independent B & Y, cavity resonance tracking,...
- Complete the move from DSAs - poor data access, network storms, free for studies,...
- LF and HF data acquisition and software to be almost identical
- Data flow for both systems same as HF Schottky in last run - proven, reliable,...
- Data will be available continuously for 6 ch of HF and 4 ch of LF Schottky/PLL

Mei Bai - Applications of RHIC AC Dipole(s)



- spin manipulation
 - spin flipping ✓
 - measure spin precession tune
- excite long-lasting coherent oscillation for
 - linear optics measurement ✓
 - linear coupling measurement ✓
 - non-linear resonance measurement
- specs

Application	BmL [Gm]	Field Orientation	Center frequency [kHz]	Tunnin range [kHz]	Duty time
Spin flipping	≥ 100	Hori	37.5	± 1.0	24s
Linear optics Measurement	78	Vert Hori	64.0	—	40ms
Non-linearity Measurement	380	Vert hori	64.0	± 1.0	≥ 80 ms




The current readback of RHIC vertical ac dipole and its spectrum. The ac dipole was configured to resonate at 36.6 kHz for the spin flipping.

Head-Tail, Button, and QMM



- # DAQ parallels present WCM system - dedicated fast scopes with deep memory, triggers from V124 or ?, local processing to spare the network
- # Head-Tail
 - Use existing 1m striplines? Install additional?
 - Useful for Instabilities, Chromaticity parasitic to ARTUS,...
 - Multiband transverse spectra
- # Buttons - Broadband transverse spectra
- # Quadrupole Moment Monitor
 - Difference of sum modes from resonant moveable/PLL BPM? New PU?
 - Resonant pickup above coherent spectrum gives needed sensitivity?
 - Injection Matching, Studies (quadrupole kicker), Schottky spectra,...

Example of Head-Tail



Triggers and Buffering



- # Want fast ‘diagnostics event’ trigger (generally on the ramp)
- # Possible triggers
 - Time domain
 - BLMs, Coherence monitor(s), Head-Tail, Buttons,...
 - Emittance growth - IPM, Schottkys,...
 - Frequency Domain - Longitudinal
 - WCM - broadband
 - LF and HF Schottky - ‘receiver’ at revolution line; narrowband
 - Frequency Domain - Transverse (Beam Transfer Functions!)
 - Buttons - broadband
 - Head - Tail - multiband
 - Schottkys - narrowband
- # For fast events necessary to buffer - a few seconds (at worst) is practical
- # What to buffer? - WCM, buttons, head-tail, QMM?, Schottky’s

Unified Data Display



- # Three monitors (ie CRTs) - two frequency domain, one time domain
- # Zoomable windows
- # Frequency domain - broadband and narrowband
 - Broadband - WCM, Buttons, Head-tail...
 - Narrowband - Artus, Schottky (10 ch?),...
 - Format - spectrogram plus most recent in conventional line presentation
- # Time Domain - Tune/Chrom/Coupling and Others
 - Tune/Chrom/Coupling - BPMs, ARTUS, PLL, Schottkys,...
 - Others - emittance, coherence, WCM/DCCT, RF BPMs, QMM, BLM?,...
- # Three Conditions (with auto-switching) - Injection, Ramp, Store
- # Data management issues - update rates, replication, display speed,...
- # Recoverable from archive in this format - unzoom time scale?