

Preinjectors Status

D. Raparia

January 19, 2016
Time Meeting

EBIS Status


- EBIS running stably Au with 6,8 or 12 pulses
- EBIS is stable with 12 pulses in 6.4 sec long super cycle
- Rotary target was put back during weekend, now running at slow rate

Linac Status


- Linac running in good condition for BLIP at 117 MeV and 152 uA with raster

- Need 1.5 hours to change 8616 tubes in mod 2

Average Current for BLIP



Beamline Location



1. plunging harp
2. 6.50'' collimator
3. bellows (3x)
4. raster magnet
5. steering magnets (existing, 2x)
6. 4.00'' collimator

7. plunging harp (opposite side)
8. laser profile monitor
9. beam current transformers (2x)
10. beam position monitor
11. beryllium window (existing)
12. aluminum bellows (existing)